

A blue and white braided suture-tendon with a curved needle, coiled on a light blue surface. The needle is positioned at the top left, and the suture-tendon extends to the right and then loops back down to form a coil in the center.


SUTUR'LINK®

Suture -Tendons  
and Ligaments Reintegration

# SUTUR'LINK®

## Storage

SUTUR'LINK® is supplied sterile, ready to use in operating rooms. It must be kept in its original unopened packaging, in a dry, clean place away from light and at room temperature.



## Product

SUTUR'LINK® is composed of a thread crimped on stainless steel needles.

It is indicated to repair or reinforce ligaments, closure and/or ligation of soft tissues, and tuberosity's reinsertions.

- Thread: USP  $\frac{3}{4}$ , L = 900 mm
- Material: 87,5 % polyethylene / 12,5 % polypropylene
- Round needle:  $\frac{3}{8}$  length 25 mm
- Triangular needle:  $\frac{1}{2}$  length 40 mm

SUTUR'LINK® is sterilized with ethylene oxide. Single use. Do not re-sterilize. For any further information, please refer to the IFU.



### DESIGNATION

Suture -Tendons and Ligaments Reintegration SUTUR'LINK®

### REFERENCE

T362034



Classe IIb



#### Head offices:

8, rue du Corps Franc-Pommiès  
65500 VIC en BIGORRE (France)  
Tél (33) 5 62 96 88 38  
Fax (33) 5 62 96 28 72

#### Administration and supply-chain offices:

Z.I de la Herray  
65500 VIC en BIGORRE (France)  
Tél (33) 5 62 96 88 38  
Fax (33) 5 62 96 28 72

#### Sales and engineering departments:

Z.I. de Montredon - 11, rue Apollo  
31240 L'UNION (France)  
Tél (33) 5 34 25 10 60  
Fax (33) 5 34 25 27 39

Distributed by

Manufactured by

